

LASER SPECTRUM ANALYZER

772B-MIR



Spectral analysis of pulsed IR and Mid-IR lasers.


The 772B-MIR Laser Spectrum Analyzer combines proven Michelson interferometer technology with fast Fourier transform analysis to characterize the spectral properties of CW and pulsed lasers that operate from 1 to 12 μm . When used with CW lasers, the model 772B-MIR operates and performs just like Bristol's 771B-MIR Laser Spectrum Analyzer. What makes the 772B-MIR system unique is that it employs a sophisticated algorithm to enable the measurement of pulsed lasers that have a repetition rate as low as 50 Hz.

When used with pulsed lasers, a single scan of the Michelson interferometer does not generate an interferogram that is sufficient to convert to a spectrum. Therefore, data from a number of laser pulses are collected from multiple scans of the interferometer to "build" a more complete interferogram.

The 772B-MIR Laser Spectrum Analyzer has a spectral resolution of 4 GHz and wavelength accuracy of ± 10 parts per million ($\pm 0.08 \text{ nm}$ at $8 \mu\text{m}$). An optical rejection ratio of about 15 - 20 dB is achieved assuming a sufficient number of laser pulses ($\geq 30,000$) are used to generate the interferogram.

Key Features:

- Spectral analysis of pulsed and CW lasers.
- Operation from 1 to 12 μm .
- Spectral resolution of 4 GHz.
- Wavelength accuracy of ± 10 parts per million.
- Continuous calibration with a built-in wavelength standard.
- Optical rejection ratio as high as 20 dB.
- Free-space aperture input with visible alignment aid.
- Straightforward operation with a PC using USB or Ethernet.
- Display software provided to control measurement parameters, display spectra, and report wavelength data.
- Five-year warranty covers all parts and labor.



Spectrum of QCL operating at a pulse rate of 50 Hz and pulse length of 50 ns

F.E.S. Co., Ltd.

1000/24, P.B. Tower, 8 floor, Sukhumvit 71 Rd., North Klongtan, Wattana, Bangkok 10110

TEL: 02-064-4050 or 02-064-4051 | FAX: 02-010-4262 | Email: info@fesupply.com

www.fesupply.com | www.thaifiberstore.com

Line Official



Facebook



SPECIFICATIONS

772B-MIR

LASER TYPE	Pulsed (repetition rate > 50 Hz, duration > 20 ns)	CW (or quasi-CW with repetition rate > 50 kHz)
WAVELENGTH		
Range ¹	1 - 12 μm	
Accuracy ^{2,3}	± 10 parts per million $\pm 0.08 \text{ nm @ } 8 \mu\text{m}$ $\pm 0.0125 \text{ cm}^{-1} \text{ @ } 1250 \text{ cm}^{-1}$ $\pm 375 \text{ MHz @ } 37.5 \text{ THz}$	± 0.75 parts per million (1-5 μm) ± 1 parts per million (5-12 μm)
Spectral Resolution ⁴	4 GHz	
Calibration	Continuous - built-in standard HeNe laser	
Display Resolution	8 digits	
Units ⁵	nm, μm , cm^{-1} , GHz, THz	
OPTICAL REJECTION RATIO ^{6,7}	10 - 20 dB (dependent on number of pulses acquired)	> 30 dB
MINIMUM INPUT POWER ⁷	0.01 - 13 μW	
MAXIMUM INPUT POWER	100 nJ (10 ns duration)	10 mW
MEASUREMENT TIME ^{7,8}	Approximately 2x time required to collect chosen number of pulses, but not less than about 10 s	< 2 s
INPUTS/OUTPUTS		
Optical Input ⁹	Collimated beam, 3 mm diameter aperture, visible tracer beam to facilitate alignment	
Instrument Interface	USB and Ethernet with Windows-based display program Library of commands (SCPI) for custom and LabVIEW programming using any PC operating system	
COMPUTER REQUIREMENTS ¹⁰	PC running Windows 10/11, 1 GB available RAM, USB 2.0 (or later) port, monitor, pointing device	
ENVIRONMENTAL ⁷		
Warm-Up Time	None	
Temperature Pressure Humidity	+15°C to +30°C (-10°C to +70°C storage) 500 - 900 mm Hg $\leq 90\%$ R.H. at +40°C (no condensation)	
DIMENSIONS AND WEIGHT		
Dimensions (H x W x D) ¹¹	7.5" x 6.5" x 15.0" (191 mm x 165 mm x 381 mm)	
Weight	14 lbs (6.3 kg)	
POWER REQUIREMENTS	90 - 264 VAC, 47 - 63 Hz, 50 VA max	
WARRANTY	5 Years (parts and labor)	

- (1) Capable of operation to 14 μm . Performance and specifications are not guaranteed beyond 12 μm .
 (2) Defined as measurement uncertainty, or maximum wavelength error, with a confidence level of $\geq 99.7\%$.
 (3) Wavelength axis is calibrated to system's accuracy specification.
 (4) Defined as the measured full width at half maximum intensity (FWHM) of an infinitely narrow optical signal.
 (5) Data in units of nm, μm , and cm^{-1} are given as vacuum values.
 (6) Acquisition of 30,000 pulses will result in an optical rejection ratio of about 15-20 dB.
 (7) Characteristic performance, but non-warranted.
 (8) For example, approximately 60 s for 30,000 pulses from a 1 kHz pulse rate laser.
 (9) Required beam height is 5.4 ± 0.25 ".
 (10) For use with Windows-based display program. Interface with SCPI can be done using any PC operating system.
 (11) Instrument height is adjustable (7.25 ± 0.25 ") for alignment purposes.



Bristol Instruments reserves the right to change the specifications as may be required to permit improvements in the design of its products. Specifications are subject to change without notice.

